Carrera®

Art. 9320231

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**Upper arm blood
pressure monitor
(battery operated)**

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


Contents

Use as prescribed	30
1. Safety instructions	32
2. Facts about blood pressure	34
Fluctuations in blood pressure	
Why is it important to monitor blood pressure regularly?	
3. The measurement process	36
4. Description of the appliance	37
5. Preparing the appliance for measurement	40
a) Inserting/replacing batteries	
b) Setting date and time	
c) Putting on the pressure cuff	
d) Important notes before measuring blood pressure	
e) Correct posture during measurement	
6. Carrying out the measurement	47
a) Switching on the appliance and automatic inflation of the cuff	
b) Saving measured values	
c) Switching off the appliance	
d) Automatic switch off	
e) Calling up the memory	
7. Description of display symbols	51
8. Specification	52
9. Cleaning and storage	53
10. Managing measured values with the PC	54
11. Disposal and environmental protection	55
12. Conformity	56
13. Guarantee and service	57

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Use as prescribed

The device is intended to be used for measuring the blood pressure on adults.

It is not intended as a substitute for an examination by a doctor. The device is intended exclusively for private use and therefore requires no subsequent calibration.

Dear Customer,

We would like to thank you for purchasing this blood pressure monitor.

It allows you to measure your blood pressure any time and anywhere precisely, quickly and easily on your upper arm.

When you use the appliance, however, you will note that there are some conditions that must be met if you wish to obtain precise, reproducible measurement results. They are explained in detail in the following.

This blood pressure monitoring device is equipped with a USB connection to save measured values on your PC.

The included software offers you the option to manage your analyzed data and view this data in various charts and graphical diagrams.

These instructions will familiarize you with your new appliance and simplify use.

Please read them thoroughly and observe the safety instructions in order to avoid accidents, injuries and damage.

Please keep these instructions in a safe place for later use and include them when transferring the appliance to a third party.

1. Safety instructions

The following words and symbols are used to indicate:

Warning!

The following safety instructions must be observed.

This symbol is used for special information about the use of the appliance and conditions of use.

- It is essential that you follow these instructions in order to obtain precise measurement results.
- This appliance does **not in any way** replace a medical examination. You should always consult your doctor if the measurement results are outside the normal range.
- With diabetics, people with pacemakers and those suffering from heart rhythm disturbances or arteriosclerosis and after the consumption of medication, blood pressures may vary from those measured by healthy patients. However, it is essential that you consult your doctor.
- Do not store the appliance in a place where it is exposed to direct sunlight, dust or moisture. You should store the blood pressure monitor and the batteries out of reach of children.

- Do **not** drop the device. It should be protected against sudden impact and vibration.
- Do **not** immerse the blood pressure monitor into water.
- Remove the batteries if you do not use the device for a month or longer.
- The appliance is intended solely for private use. Therefore, it is not covered by any calibration requirement.
- It is not suitable to measure childrens' blood pressure by means of this device with the provided cuff.
- Switch the appliance on only if the cuff is correctly put around the upper arm.
- Contact our Service Center (see page 57) in case of defects.
- **This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities or lack of experience and knowledge, unless they have been given supervision or instruction concerning the use of the appliance by a person responsible for their safety.**
Children should be supervised to ensure that they do not play with the appliance.

2. Facts about blood pressure

Blood pressure is the force which moves the flow of blood around the body. This also exerts pressure on the walls of the arteries and this pressure can be measured.

The arterial blood pressure depends on the heart's pumping cycle.

The highest pressure in this cycle is known as **systolic** blood pressure, and the lowest as **diastolic** blood pressure.

Your doctor needs both the **systolic** and the **diastolic** blood pressure to be able to assess your circulation.

Fluctuations in blood pressure

Individual blood pressure fluctuates depending on the time of day and season.

It is influenced by various factors including physical activity, fear, nervousness, etc. These fluctuations are even more marked in people with high blood pressure. Blood pressure is usually lower during sleep and before getting up and then rises.

Why is it important to monitor the blood pressure regularly?

The monitoring of blood pressure by a doctor as part of a medical examination can cause the person being examined to feel nervous and thus cause a rise in blood pressure. Our blood pressure also depends on various other factors and can therefore not be assessed on the basis of a single measurement.

Even in the short term, it can fluctuate considerably which is why repeated measurements can provide different results. If the patient records his blood pressure over the course of the day, it can therefore provide valuable support for the doctor.

Blood pressure measured directly in the morning after getting up and before breakfast when the person is at rest is described as the fundamental blood pressure. In practice, this is very difficult to measure.

In order to record it as precisely as possible in this environment, however, it makes sense to carry out the monitoring at home.

The following new standard values for blood pressure have been defined by the World Health Organisation (WHO) without reference to age:

Assessment	Systolic (mm Hg)	Diastolic (mm Hg)
Optimal blood pressure	< 120	< 80
Normal blood pressure	120–129	80–84
High-normal blood pressure	130–139	85–89
Mild hypertension (level 1)	140–159	90–99
Medium hypertension (level 2)	160–179	100–109
Severe hypertension (level 3)	> 180	> 110

Source: WHO, 1999

3. The measurement process

This compact, fully automatic upper arm blood pressure monitor measures blood pressure using the oscillometric method.

The appliance has an upper arm cuff which can be put on easily. At the press of a button it measures your blood pressure and your pulse and shows the values measured on an easy-to-read digital display.


The appliance can save 100 measurements for one user in total. It can also be set up for up to 2 users, in which case 50 measurements can be saved for each one.

Therefore, the appliance is suitable for a quick and easy blood pressure measurement at home.


4. Description of the appliance

Overview of operating elements

Measurement device

Display symbols

Explanation of display symbols

Symbol	Range	Systolic	Diastolic
	Hypertension	> 140 mmHg	> 90 mmHg
	Borderline	130-139 mmHg	85-89 mmHg
	Normal	< 129 mmHg	< 84 mmHg
	Hypotension	< 90 mmHg	< 50 mmHg

Accessories

- ① Blood pressure monitor
- ② Blood pressure cuff with air tube
- ③ Storage bag
- ④ CD with software and user manual
- ⑤ USB cable
- ⑥ 4x LR6 (AA) 1.5 V batteries

5. Preparing the appliance for measurement

a) Inserting/replacing batteries

1. Put the appliance on its front and pull off the battery compartment lid.

Now, insert the batteries according to the position shown in the compartment.

Make sure that the battery poles (+/-) are the right way round.

2. In case of missing display and/or display of empty batteries symbol all batteries must be replaced.

 Symbol indicates empty batteries

i NOTE: Insert the new batteries within 30 seconds after removing the old batteries, since otherwise all the saved data will be deleted.

i Important: In order for the appliance to work properly, the battery compartment cover must be firmly in place. Remove the batteries from the appliance if you do not use it for any length of time. Batteries should be disposed of according to the regulations (see point 11).

b) Setting date and time

The default value of this function is *1/01 12:00*

1. To change the date and time, press the Adjust **T**-button for about 4 seconds. The "Year" starts to blink on the LCD.


2. If you press Time+ **M1** or Time– **M2** button you increase or decrease the number.
3. The flashing figure moves on one digit when the Adjust **T** button is pushed. The flashing display changes from month to day, hour and minute, if you press the time button repeatedly. If you press Time+ **M1** or Time– **M2** button you increase or decrease the number.
4. End the setting by pressing the Adjust **T** button after you have changed the date.
The number indicating the year disappears behind the time and is displayed for approx. 1 second when briefly pressing the Adjust **T** button.
5. To cancel the setting, leave the appliance without pressing any buttons. Ten seconds later the display will show the original display and the current time.

c) Putting on the pressure cuff

1. Insert the air hose plug on the pressure cuff firmly into the left socket on the measuring device.
2. Open the cuff out into a ring and pull it over the left wrist up to the left upper arm (without clothing).
Position the pressure cuff on the upper arm, around 2-3 cm above the elbow, as shown in the drawing on the cuff.
3. Pull the cuff end, which is taken through the metal clip, outwards across the upper arm and fix it with the Velcro closure.

i NOTE:

- Put the pressure cuff on your **left** upper arm in order to obtain precise measurement results.
- Only use the original cuff, since other cuffs may adversely affect the measurement.
The hose connection must face the body and should be in the middle of the bend of the arm.
- Do not pull the pressure cuff too tightly.
- The cuff clip should not be over the artery, as this could cause incorrect measurements.

- The smooth side of the cuff should face inwards.
- Do not pull on the pressure cuff hose.

Warning!

With this appliance, the compressed air is automatically expelled quickly once the cuff pressure has reached or passed 320 mmHg.

To prevent the cuff tightening too much around the upper arm in the event of a faulty appliance, the user must remove the cuff from the main unit immediately.

i d) Important notes before measuring blood pressure

1. Carry out the measurement somewhere quiet and in a relaxed position.
2. Stress raises the blood pressure! You should therefore not measure your blood pressure in stressful situations.
3. In rare cases, errors may occur preventing correct measurement in people with an extremely weak or irregular pulse.
4. It is normal for blood pressure to change in the course of the day. It is also influenced by various factors such as smoking, alcohol consumption, the use of medicines and physical activity.
5. Pacemakers and changes in the arm veins can cause the measurement to be imprecise.
6. The blood pressure measurement can also be influenced by the position of the person examined, his/her physical condition and other factors.
7. Stay calm during the measurement, and do not talk since this could cause an incorrect measurement.
8. Wait five minutes before repeating the measurement so that your circulation can return to normal. Depending on your physical characteristics, a longer or shorter interval may be required between measurements.

9. Blood pressure measurements must be interpreted by doctors or trained medical staff who are familiar with your medical history. By using the monitor regularly and recording the results for your doctor, you can inform him about ongoing tendencies shown by your blood pressure.
10. To ensure comparability, always carry out the measurement at the same time of day, if possible, since blood pressure can fluctuate considerably.
11. Please keep this appliance away from high electromagnetic fields such as TV and radio appliances, etc.

e) Correct posture during measurement

The measurement values depend on your physical posture.

- You should only measure your blood pressure when you are sitting down.
Always make sure that you are sitting in a comfortable, relaxed position.
- Place your elbow on a table or other object.
If necessary, use a suitable base to position your arm in such a way that the pressure cuff is at the same level as your heart.
Your blood pressure will vary depending on the level of the cuff in relation to your heart.
You must stay in this position until the measurement is completely finished.
- Slightly open your hand and relax.

i **NOTE:**
Your heart is slightly below your armpit.

Important if you want to obtain correct measurement results at all times:

During the measurement:

- Do not move, and do not exert any pressure on your left wrist or left hand.
- Do not talk, sneeze or cough.
- Do not carry out any measurements when you are standing up.

If you are repeating the measurement, wait five minutes between measurements.

In this time, the congested blood vessels can return to their normal condition before the next measurement is carried out.

The measuring results may not be reproducible at all times when deviating from the steps described here!


6. Carrying out the measurement

Before carrying out the measurement, please read the following points of Chapter 5 through completely:

- c) Putting on the pressure cuff
- d) Important notes before measuring the blood pressure
- e) Correct posture during measurement


a) Switching on the appliance and automatic inflation of the cuff

- Press the **Start** button of the appliance.
- When the appliance is switched on, all the display symbols appear for approx. one second.
- The display symbols disappear again, during display of systolic blood pressure a blinking \square appears.
The pressure cuff is pumped up within a few seconds.

- You will be able to see the pump pressure raise in the display. After inflation has been completed and a short standstill, the pressure drops gradually.
- After inflation, the ♥ symbol appears to indicate that measurement is ongoing. Inflation of the cuff again is perfectly harmless and normal for some users.

- The reduction of pressure starts by the solenoid valve opening briefly. Once the pressure has been slightly reduced and pulse recognised, the pressure reduction is synchronised with the pulse of the person examined and the pressure is released completely.
- If the **Start** button is pressed during the build-up of pressure, the appliance is switched off and the air expelled from the cuff.
- When the measurement is completed, the pressure release is taking place. The systolic pressure appears in the upper display area, the diastolic pressure in the middle area, and the pulse in the lower display area.

These measured values are displayed as long as the **Start** button is pressed again.

They get lost if they are not saved.


After 3 minutes the values are automatically deleted and the appliance changes to the "Time" display mode.

If you repeat the measurement, wait 5 minutes between measurements.

b) Saving measured values

Should the measured values be saved, it is necessary to press the **M1** or **M2** button.

Then, the values are assigned to the respective user and lodged in his memory area.

When pressing the **M1** (or **M2**) button the **M/R** sign flashes and the measured value is assigned to the data memory of user 1 (or 2).

If a new measurement was effected, the new measured value is lodged into the save location 1 and the value previously saved automatically moved to save location 2 etc.

The measured value of the last save location gets lost if the save capacity of 50 values per user is passed.

c) Switching off the appliance

Switch off the appliance using the **Start** button.

Date and time indication remains unchanged.

d) Automatic switch off

If the appliance is not switched off after the measurement, an automatic shutoff function switches it off after about 3 minutes.

e) Calling up the memory

- To call up an earlier measurement, press the **Time M1** or **M2** button.

The last measured value is displayed for 4 seconds. The save location number appears in symbol "Book".

Then, display automatically returns to the previous time indication if no entry is made.

- The save number increases by one level if the **Save** button is pressed again whereas the appropriate measured value is displayed.

7. Description of display symbols

Display symbol	Status/cause	Remedy
	Appears during measurement and flashes once the pulse is detected.	Ongoing measurement – stay calm!
	Display when batteries are empty.	Replace all 4 batteries with new ones.
$E1$ -1	Pressure raise is defective. Cuff is not correctly put on or leaky and/or appliance is broken.	Check cuff and fixing. Repeat measurement after 5 minutes.
$E1$ -2	Measurement was incorrect. User has moved, is not calm or has the wrong position.	Stay calm! Repeat measurement after 5 minutes.
$E1$ -3	Pressure decreases in cuff too fast. Cuff is loose or device is defective.	Put cuff on again. Repeat measurement after 5 minutes.
$E1$ -4	Abnormal measured values or appliance is defective.	Maybe, user has problems with health and should consult the doctor.
$H 1$	Overpressure in measurement system. Tube blocked or device is defective.	Check tube and cuff. Repeat measurement after 5 minutes. Consult service.

8. Specification

Specification and model number	Fully automatic upper arm blood pressure monitor 9320231 (BP-600U)
Protection against electric shock	Constructed according to type BF 
Display system	Digital LCD display system SYS, DIA, PULSE, TIME, DATE, USER
Measurement method	Oscillometric method
Power supply	LR6 (AA) 4x 1.5 V batteries
Measurement range	20–280 mmHg (pressure) 40–200 beats/minute (pulse rate)
Precision	± 3 mmHg (pressure) ± 5 % (pulse)
Initial pressure setting	180 mmHg
Pressure system	Automatic pressure build-up with micropump
Fast pressure release	Solenoid valve
Air expulsion	Automatic linear air expulsion with solenoid valve
Memory	100 measurements for one user 50 measurements each for two users
Automatic shutoff	Automatic shutoff after 3 minutes idleness to protect batteries

Automatic build-up of pressure again	Automatic build-up of pressure again if cuff pressure is insufficient
Working life of battery	approx. 200 measurements
Operating environment	+10 °C to +40 °C, 30 %–85 % rel. hum.
Storage environment	– 5 °C to +50 °C, 15 %–90 % rel. hum.
Arm span	approx. 230–330 mm
Accessories	Storage bag, instruction manual, 4 batteries LR6 (AA), 1.5 V, USB cable, CD

Special Accessories:

Extra large cuff (circumference of approx. 33 to 43 cm).

Please order from the Service Center (see p. 57). **Art. 9320331**

9. Cleaning and storage

1. The appliance and cuff should only be wiped with a damp cloth.
The cuff may not be washed or dry cleaned.
2. Do not store the appliance where it is exposed to direct sunlight, dust or severe moisture.
3. Keep the appliance and the batteries out of reach of children.
4. Remove the batteries if you do not use the appliance for any length of time.

10. Managing measured values with the PC

1. The included CD features the user manual of the manufacturer (English) in the form of a PDF file to be viewed with your PC. Open the file with Acrobat Reader and print out the manual.
2. The software can be installed with the following Windows operating systems:
Win98SE/ME/Win2K/XP/2003/Vista32/Win7.
3. Insert the CD and click "Install Software" after the installation dialog opens (see Fig.). Follow the steps listed on screen to install the software.
4. Once the software has been installed successfully on your PC, you are able to transfer the saved measured values from the blood pressure monitor to the PC using the USB cable and then manage your data as desired (**PDF instructions, page 14**). Please note the respective model designation.

Note: Direct measurements on the monitor of the computer are not possible. Any data you wish to view on your PC must be transferred via USB cable first.

11. Disposal and environmental protection

Please do not throw the device into the refuse bin at the end of its service life. Enquire at your municipal offices about methods of disposing of it as electronic waste properly in an environmentally friendly manner.

Under the law, you are not allowed to throw batteries into the domestic waste. Please dispose of batteries via the battery return system or using one of the special collector containers provided in shops.

Packaging materials are resources and can be reused. Dispose of these in an environmentally friendly manner and as stipulated by local rules and regulations.

HINWEIS:

Within the scope of product improvements we reserve the right to make technical and visual changes to the **Health Control Professional** and accessories.

12. Conformity

The appliance is an application component, type BF, and complies with European Directives:

93/42/EC Directive on Medical Equipment

The CE marking has been affixed to the product.

The associated conformity declaration is held by the EU representative.

INFORMATION ON THE GUARANTEE:

Please note that incorrect or improper handling of the device, failure to comply with the safety precautions that apply to the device, use of force or alterations or repairs that are carried out by a party other than the service address that we have authorized void the guarantee.

Furthermore, components that are subject to normal wear and tear or use are excluded from the guarantee.

13. Guarantee and Service

Under our terms of warranty, we guarantee that our technical appliances will provide faultless quality and function.

The guarantee period is 36 months from the date of purchase. Please keep your receipt as proof of purchase. The trader's legal guarantee obligation is not limited by this guarantee. During the guarantee period, you can send defective appliances to the Service Center (see address below), **postage paid**. Then you will receive a new appliance or get the repaired one back free of charge. **Parcels not prepaid cannot be accepted.**

Please use the original packaging in order to avoid any damages in transit and enclose a copy of the receipt. A detailed description of the complaint will help reduce repair time.

After expiry of the guarantee period, you can still send a defective appliance to the address below, also **postage paid**. Then, however, repairs made will be at your expense.

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ELRAS GmbH – Service Center –
Produced by: YA HORNG ELECTRONIC CO., LTD.